



Order Form: Barwon Lift Chair

Prepared by

Client Name

Order Date

Ref no.

OT's Name

OT's Email

Funding Option

☐

Aged Care

☐

NDIS

☐

Self

☐

Other (Specify)



CONFIGURATION

Essentials

☐

Dual Motor

☐

180kg Dual Motor ^{RC}

☐

Single Motor

☐

Manual Recline ^{RC}

☐

Comfort Pack

☐

Heat & Massage ^{CP}

☐

24 Volt Battery

☐

Tilt in Space ^{HM}

Wheels

☐

50mm Castors
(Standard)

☐

75 Mobile Castors
(+25mm Seat height)

☐

Glides

☐

Remove Castors
(in Factory)

Memory Foam

☐

Seat

☐

Legrest

☐

Backrest

☐

Armrest

☐

Hygroflex

Back Cushion Options

☐

Gluteal Recess

☐

Concave Back

☐

Blow-fill Lumbar ^{CP}

Cushion Insert

☐

None
(Standard)

☐

18x16
Low Profile

☐

18x16
High Profile

☐

16x16
Low Profile

☐

16x16
High Profile

☐

Include
Infill Plug

Controller Option

☐

Right (Standard)

☐

Left

☐

Double Pocket

☐

Controller Loop
Left / Right

Back Angle

☐

Standard

☐

Vertical

☐

Over-Vertical

Other Options

☐

Waterfall Arms

☐

Scotch Guard

☐

Flush Front

☐

Meal Tray Left/Right (default show wood removed)

☐

Leg rest Side Walls

☐

No arm
chair cover

☐

Full
chair cover

☐

Velcro armrest/
headrest covers

☐

Loose armrest/
headrest covers

☐

Fitted armrest/
headrest covers

☐

Removable Headrest **A** **B** **C** **D** **E**
(please circle one)

COLOUR AND FINISH

Fabric

Inner

Outer

Show wood

☐

Walnut

☐

Mahogany

☐

Baltic

☐

Teak

☐

Remove Show Wood



Order Form: Barwon Lift Chair

Circle required dimensions and configurations

DIMENSIONS

Size	Seat Depth (mm)	Seat Height (mm)	Back Height (mm)	Footrest Length (mm)	Scapular Height (mm)	Seat Width (mm)	Arm Height (mm)	Seat Foam (mm)
------	-----------------	------------------	------------------	----------------------	----------------------	-----------------	-----------------	----------------

Decrease Size	-100 CP		-150 CP		-36 MM		-120 CP		-150 CP		-180 CP		-50			-150		-60																						
	-25		-50		-75		-25 RC		-30 CP		-60 CP		-90 CP		-25			-50		-100		-30		-20																
Mini	455					456					636					N/A					429					540					199					N/A				
Increase Size	+25		+50		+20		+30		+60		+90		+25			+50		+30		+20																				
	+75										+120		+150 DM		+50					+60																				

Decrease Size	-100		-150 CP			-120 CP		-150 CP		-180 CP		-50			-150		-60																		
	-25		-50		-75		-25 RC		-30		-60 CP		-90 CP		-25			-50		-100		-30		-20											
A				490				476				647				N/A				440				540				180				N/A			
Increase Size	+25		+50			+20		+30		+60		+90		+25			+50		+30		+20														
	+75												+120 DM		+150 DM								+60												

Decrease Size	-100		-150 CP			-120 CP		-150 CP		-180 CP		-50			-150		-60																		
	-25		-50		-75				-25 RC		-30		-60		-90 CP		-25			-50		-100		-30		-20									
B				520				496				688				N/A				481				540				190				N/A			
Increase Size	+25								+25				+30		+60 DM		+90 DM		+25								+50				+30		+20		
	+50												+120 DM		+150 DM														+60						

Decrease Size	-100		-150			-120 _{CP}		-150 _{CP}		-180 _{CP}		-50			-150		-60			
	-25	-50	-75		-25 _{RC}		-30	-60	-90			-25			-50	-100	-30	-20		
C																				
545																				
546																				
718																				
N/A																				
511																				
540																				
190																				
N/A																				
Increase Size	+25				+25				+30 _{DM}		+60 _{DM}		+25		+50		+30		+20	
					+50				+90 _{DM}		+120 _{DM}		+50				+60			

ADDITIONAL NOTES

NOTE:

- * Reduced seat foam limits Safe working load to 90kg and reduces seat foam warranty to 2 years.
- * Seat width decrease will automatically come with waterfall arms.
- * Castors can only be removed in factory.
- * Fusion fabrics only available in 2-Tones.
- * Full chair cover and Fitted arm covers not compatible with Waterfall arms.
- * Blow-fill lumbar not available with memory foam back.

* MM - Micro Mini

- * Only the following selections are available on Micro Mini: Cushion insert, Controller option, memory foam.
- * Gluteal recess and concave back are the only available back cushion options.

* **Prescriber / OT** = not compulsory for completion

* **RC** = Remove Castors -25mm Seat Height


* **CP** = Not compatible with Comfort Pack

* **HM** = Not compatible with Heat and Massage

* **DM** = Dual Motor Only

Chair Sizing

OUR RECOMMENDED STARTING POINT



C	5'11+	180cm+
B	5'8 to 5'11	172 to 180cm
A	5'2 to 5'8	157 to 172cm
Mini	4'10 to 5'2	147 to 157cm

ENSURING CORRECT MODIFICATIONS



- Heels up against the footrest
- Light pressure on the calf
- Feet planted on the ground, avoid legs angled too far forward
- Mention gluteal recess if user unable to sit back into the chair



- Nape of the neck sitting in centre of top cushion
- Shoulder blades sitting in gap of the top roll
- Adjustments can be made manually by repositioning backrest



- Knees in line with hips
- Feet planted on the ground
- Even pressure from their backside to the back of their legs
- Use blocks to determine adjustment if feet don't reach the ground



- Elbow to wrist should be supported
- If shoulders are hunched up, arms too high
- Gap between armrest and elbow, arms too low



- Fit a fist beside each thigh comfortably
- Waterfall arms can also decrease seat width



- Heel clear of footrest to avoid pressure sores
- Ankles to be supported